



INVAMED



SHOPRIDER®

Traveso

Model S8889SXL SBN/CAN | Class 3 Scooter

Traveso

Model S8889SLSBN/CAN | Class 3 Scooter



Luxurious rustproof cab.



Two lockable doors.



Windscreen wiper.



Headlights.



Non-marking pneumatic tyres.



Wing mirror.



York Road
Bridgend Industrial Estate
Bridgend
CF31 3TB

Tel: 01656 674488
Fax: 01656 674499
Email: sales@invamed.co.uk

www.invamed.co.uk

 We are social. Follow us on Facebook
@invamedonline

*Lower limb amputees should refer to the user manual before use.
** Passenger weight, ground conditions and battery condition/size can affect distance of travel.
*** Measurement of underside of chassis.

Roma Medical reserve the right to change specifications without prior notice.



Special Features

- Luxurious rustproof cab.
- Two lockable doors*.
- Headlights.
- Rear lights.
- Indicators.
- Wide profile wheels.
- Non-marking pneumatic tyres.
- High/low speed switch.
- Captain Seat - Adjustable.
- Headrest.
- Full suspension.
- Delta tiller.
- Horn.
- Windscreen wiper.
- Windscreen washer.
- Window vents.
- Wing mirrors.

Specifications

Model No.*:	S8889XLSBN/CAN
Class (Invalid Carriage):	3
Overall Length:	1600mm / 63"
Overall Width:	730mm / 29"
Overall Height:	1620mm / 64"
Seat Width:	500mm / 20"
Safe Working Slope:	8° / 14%
Maximum User Weight:	160kg / 25st
Battery Size:	2 x 12v (36Ah - 75Ah)
Range**:	16 -32 Km / 10 - 20 Miles
Maximum Speed:	6 & 12 Kmph / 4 & 8 Mph
Wheel Size (Front):	330mm / 13"
Wheel Size (Back):	330mm / 13"
Turning Radius:	1780mm / 70"
Battery Charger:	8 amp off board
Ground Clearance***:	60mm / 2.5"
Colours:	Silver